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9000 St. Gallen / Switzerland

Tungsten Carbides Tungsten carbide burs/finishers for Dental Surgeries

Recommendations for Use and Safe Operation in Dental Practices

1. General

To prevent harm to the residual tooth structure and the tissues and teeth contiguous to the site, provide for sufficient cooling, technically, clinically and hygienically optimum instruments and materials as well as optimum instrument speed and only exert the necessary amount of pressure. The service-life of tungsten carbide instruments depends on compliance with the following recommendations for use and safe operation.

2. Recommended hygiene procedures

The tungsten carbide burs are selected (shape, size, type) according to the type of preparation to be carried out.

The ergonomic principles must be taken into account during the working procedures.

The handpieces (high-speed handpieces, contra-angles) for tungsten carbide instruments must be in perfect working conditions.

The instruments must be fully inserted and gripped in the high-speed handpiece or contra-angle with the utmost of care.

Select the correct speed to ensure optimum performance and service-life. Operate within the recommended speed ranges.

The hints covering maximum permissible speeds must be heeded.

Before placing the tungsten carbide bur/finisher in contact with the site being prepared, check that the bur is rotating at the optimum speed and that the cooling system is functioning properly. Once the preparation has been completed, remove the rotary instrument from the site and allow it to come to a standstill.

Prepare carefully and without exerting pressure. Select the pressure (0.3 – 2 N) to ensure that the speed is not reduced noticeably. The pressure exerted on the rotary instrument is a decisive factor for the working procedure and the successful outcome of the preparation. Under no circumstances should the burs be allowed to jam.

To prevent harm to the tooth structure, pulp and adjacent fillings, cool by spraying with sufficient amounts of water. We recommend high-speed handpieces with three jets which spray the coolant along the entire length of the instrument's cutting surface. Preparing a tooth without using sufficient coolant – minimum: 50ml/min. and incorrect application of the coolant (spray diverted, inadequate ejection) may adversely affect the final result.

Only use impeccable tungsten carbide instruments. Instruments which are blunt, damaged, bent or no longer concentric, must be segregated immediately and disposed of carefully.

Whenever possible, the tungsten carbide instruments should be stored in their original packaging to enable them to be identified and traced.

Recommended Hygiene Procedures for Dental Practices

1. General

Protection against infection is especially important to patients and surgery staff. Disinfection and/or sterilization of the equipment and instruments used are an absolutely essential part of a safe preparation technique.

- Tungsten carbide instruments must be sterilized prior to use.
- The following hygiene procedures are recommended before and after use.

Date of issue:

06-2018

Only supplied to dentists and dental laboratories or upon their instructions.

2. Recommended hygiene procedures

Cleaning, disinfecting, sterilizing and storing tungsten carbide instruments:

Valid methods for sterilizing and cleaning burs are:

- *Autoclaving:* Dip instruments in a corrosion inhibitor such as 1% sodium nitrite, prior to the sterilization cycle. Operate at full cycle with a dwell at 132°C minimum for 12 minutes.
- *Dry Heat:* Expose the instruments to a temperature of 190°C for 6 minutes.
- *Cleaning with Ultrasound:* Instruments may be ultrasonically cleaned by inserting them in instrument holders (or blocks) to prevent damage by rubbing against other surfaces or each other. A cycle of 5 minutes is recommended, using a general purpose cleaner.

Caution:

- Use sterilizing devices according to the manufacturers recommended procedure. It is the responsibility of the user to ensure that sterilization is effective.
- Avoid cold disinfectant solutions. These solutions do not sterilize and may contain corrosive agents that degrade the performance and strength of the carbide instrument.
- If necessary, clean by mechanical means in advance (brush the instrument blades with a brush with nylon or brass bristles). Check for any remains and, if necessary, clean again.
- Dental instruments made of carbide and steel should be dried and stored in a moisture free environment. This avoids the possibility of corrosion which could weaken the joint of carbide and steel.

Warning:

This product may contain nickel and should not be used on persons with known allergic reactions to nickel or nickel alloys.

Please note:

To prevent spreading infectious germs, instruments contaminated with blood or saliva should be disposed of carefully.

Recommended Speeds

Head diameter in 1/10 mm	Speed range (r.p.m.)
005 – 006	53.000 – 105.000
007 – 008	40.000 – 80.000
009 – 010	30.000 – 60.000
012 – 014	24.000 – 48.000
016 – 018	19.000 – 37.000
021 – 023	15.000 – 30.000
025 – 027	13.000 – 26.000
029 – 031	12.000 – 24.000

Marking

The **LOT** number is shown on the package.

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